

Deluxe Folding Walker, Two Button with 5" Wheels

HCPCS:



FIGURE A



FIGURE B




FIGURE C

- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Sturdy 1" diameter aluminum construction ensures maximum strength while remaining lightweight
- Newly designed rear glide cap allows for walker to slide easily and smoothly over most surfaces (Figure A)
- Vinyl-contoured hand grip
- Includes wheels (Replacement wheels item # 10109)
- Polypropylene (PP) and PVC wheels allow walker to roll easily and smoothly over irregular surfaces
- The 10210-4ASM Walker comes fully assembled out of the box
- Piggy back label is supplied on product hangtag

PRODUCTS

Description

10210-1	Walker Adlt 2but w/5 whl.1/ea
10210-4ASM	Walker,2BtnAdult,5"WhlsAttached,4/cs
10211-4ASM	Dlx.Walker,2BtnJunior,5"WhlAttached,4/cs

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Handle Height	32"-39" 25"-32.25"
Weight Capacity	350 lb
Carton Shipping Weight	30 29 9.5 lb
Casters	5"
Overall Width	24"
Weight	7.5 7.5 7 lb