

Lightweight Bariatric Homecare Bed

HCPCS: 15300LW: E0301|15302LW and 15303LW: E0302,E0302



Drive is leading the way in Homecare beds by introducing the next generation of Bariatric support. Introducing Drive DeVilbiss' new Lightweight Bariatric Beds. With frames up to 100 lbs lighter than most current Bariatric offerings, these beds are easy to maneuver, transport and assemble. Without compromising strength, these durable lightweight frames have been tested and verified to maintain the same safe working load capacity of up to 1000 lbs! (Weight capacity varies by size). Available in 42", 48" and 54" widths gives patients the space they need to recover and providers the ability to deliver and assemble with ease.

- Lightweight frame design offers providers easier and more manageable bed deliveries and assemblies
- Allows for additional space for patient comfort and recovery as well as caregiver efficiency
- Product Safety & Performance Standards: Meets IEC 60601-2-52; 60601-1; 60601-1-2; 60601-1-11
- Safe working load capacity of 600 lbs
- Full electric bariatric homecare bed with dual motor assembly for head and frame adjustability
- Ultra lightweight frame is easily set up by 1 individual; no links to connect and no tools required
- Square tube steel frame slat deck provides durable strength, no-sag mattress support, and is easy to clean
- Scratch resistant wood grain panel bed ends include hi lo motors for frame height adjustability
- UL approved motor is self contained for reduced weight and noise and can be easily installed and removed with patient in bed
- Reduce service calls with large glow in the dark hand pendant that is easy-to-use for adjustments; LED indicators available for power and lockout included
- Battery Back up accessory available (#15302LWB) and sold separately
- Half Rails accessory available (15021BV-LW) and sold separately
- Rail accessories meet FDA and ASTM entrapment guidelines and standards

PRODUCTS	Description
15300LW	Lightweight Full Elec 42" Bariatric Bed
15300LW-HR	Lightweight 42" Bari Bed w/Half Rails
15302LW	Lightweight Full Elec 48" Bariatric Bed
15302LW-HR	Lightweight 48" Bari Bed w/Half Rails
15303LW	Lightweight Full Elec 54" Bariatric Bed
15303LW-HR	Lightweight 54" Bari Bed w/Half Rails

SPECIFICATIONS	Description
Warranty	5 years on Motor and Hand Control Limited Lifetime on Welds and Frame 1 year on all other parts and components
Battery	Backup available but not included
Casters	5" casters (2 locking and 2 unlocking)
Color	Brown

ACCESSORIES**Description**

15302LWB	LW Bari Bed Accessory battery pk,48/54, BX
15021BV-LW	Half Side Rail,42in/48/54 LW Baribed,1 pr, BX