

Bariatric Aluminum Two Button Folding Walker with Wheels

HCPCS:




FIGURE A



- For Bariatric use
- Steel legs and side braces
- Precision design provides additional strength while adding minimal weight
- Wider and deeper frame design when compared to standard walkers
- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability
- Vinyl-contoured hand grip

PRODUCTS	Description
10220-1WW	Walker, AdltHD, 2Btn w/Whls, 1/cs

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	500 lb
Height	32" - 39"
Carton Shipping Weight	24 lb
Overall Width	27"
Weight	11 lb