

Silver Sport Recliner with Vinyl Upholstery

HCPCS: K0001 + E1226 + K0195,K0001 + E1226 + K0195



- Vinyl upholstery is comfortable, lightweight and easy-to-clean
- Hydraulic reclining mechanism allows for infinite adjustments up to 180°
- Arms include support extension for reclining position
- Carbon Steel frame with silver vein finish is attractive, and easy-to-maintain
- Composite Mag-style wheels are lightweight and maintenance free
- 8" front casters
- Padded armrests provide added patient comfort
- Wheels are set back on frame to prevent tipping
- Precision sealed wheel bearings in front and rear ensure long-lasting performance and reliability
- Rear anti-tippers come standard
- Headrest extension with cushioned head immobilizer come standard
- Includes swing-away elevating leg rests with Metal foot plates
- Comes with push to lock wheel locks

PRODUCTS	Description
SSP16RBDDAV	SilverSport16"ReclBack,DetDesk,Vnl,1/ea
SSP16RBDFAV	SilverSport16"ReclBack,FullArm,Vnl,1/ea
SSP18RBDDAV	SilverSport18"ReclBack,DetDesk,Vnl,1/ea
SSP18RBDFAV	SilverSport18"ReclBack,FullArm,Vnl,1/ea
SSP20RBDDAV	SilverSport20"ReclBack,DetDesk,Vnl,1/ea
SSP20RBDFAV	SilverSport20"ReclBack,FullArm,Vnl,1/ea

SPECIFICATIONS	Description
Warranty	Limited Lifetime
Product Weight	59 lb 60 lb 61 lb
Product Dimensions	50" (L) w/ riggings x 23" (W) x 50" (H)
Material	Steel frame Vinyl upholstery
Seat Dimensions	16" (W) x 16" (D)
Weight Capacity	16" - 18" 300 lb 20" 350 lb
Seat To Floor Height	20.5"
Carton Shipping Weight	66 lb 67 lb 67.5 lb
Folded Dimensions	50" (L) x 12" (W) x 50" (H)
CA Prop 65	Yes
Back Height	21" 31" w/ headrest
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]

ACCESSORIES	Description
STDS856	Seat Belt,60"Velcro-type,1/ea, EA
STDS850	Seat Belt, Clasp Type,1/ea, EA
STDS855	Seat Belt,60" Auto-type,1/ea, EA
STDS1040S	Cup Holder, Swivel,1BX,6/CS, BX
STDS806	Anti-Fold/Theft Device for w/c,1/ea, BX
STDS818	Anti-Tippers,rear,SentraRecliner,1pr/bx, BX
STDELR-AL	ArticulatingElev.Legrest, Std,1pr/cs, CS
LK3JELR	ELR Set,Black,1Pair eq EA, PR
WASR	Limb Support, Right,1/cs, CS