

Premium Guard Gel Overlay

HCPCS: E0185



FIGURE A

- Best in class gel overlay optimizes pressure redistribution and is now available in up to 500 lbs weight capacity
- 15 gel bladders cover entire patient contact area for the prevention and treatment of pressure ulcers
- Full length "head-to-toe" gel bladders minimize "bottoming out" and "pitting"
- New convoluted foam topper vastly improves patient comfort and air circulation
- 4 elastic corner straps secure overlay to mattress
- Urethane coated nylon cover is water resistant and vapor permeable
- Vinyl base is stain resistant and waterproof
- Bladders contain an aqueous non-toxic gel and polystyrene beads to minimize gel migration
- Meets flammability standard TB117

PRODUCTS	Description
14893	PremGuardGel/FoamOverlay34x76x3.5in,1/cs
14901	Prem.GuardGel/FoamOverlay42x76x3.5",1/cs
SPECIFICATIONS	Description
Warranty	Foam: 18 Months Cover: 18 Months
Weight Capacity	300-500 lbs
HCPCS	E0185
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]
ACCESSORIES	Description
14893C	Prem.Guard Gel,3"Foam Overlay,Cover,1/ea, EA