

3-Position Recliner

HCPCS:



FIGURE A



FIGURE B



FIGURE C

- 3-positions: upright, deep recline and elevated leg rest/footrest
- Retractable lock bar secures the chair in the desired position
- Self-adjusting suspension in the seat prevents bottoming out
- Side panels "pop off" for easy cleaning and maintenance
- Large, blow-molded tray locks in 5 positions and can be stored on the side of the chair
- 5" casters (2 with locks) make it easy-to-move and secure the recliner
- Wide, comfortable armrests are angled to provide support in any reclined position
- Neutral-colored, sturdy steel frame is easy-to-clean and maintain
- Moisture barrier on seat prevents seepage
- Comfortable built-in headrest
- Foam and vinyl meet California Technical Bulletin 117
- Standard

PRODUCTS	Description
D574-CHAR	3 Position Standard Charcoal
D574-CHOC	3 Position Standard Chocolate
D574-TAN	3 Position Standard Tan
SPECIFICATIONS	Description
Warranty	Limited Lifetime on Frame, 1 Year on Upholstery and Moving Parts
Seat Dimensions	21" (H) x 19" (W) x 19" (D)
Weight Capacity	250 lbs
Carton Shipping Weight	110 lbs.
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]