

Side Walker

HCPCS: E0135



FIGURE A



- Designed for individuals with the use of only one hand or arm
- Lighter than a walker and more stable than a cane
- Folds easily with one hand (Figure A)
- Ideal for users who do not require a walker but need a wider base for support

PRODUCTS	Description
10240-1	SideStyle Walker,1/ea
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Weight Capacity	300 lb
Height	29.5" - 36.5"
Carton Shipping Weight	5 lb
Weight	4 lb