

Hydraulic, Deluxe Silver Vein Patient Lift

HCPCS: E0630



- Silver Vein steel construction
- 6 point Swivel Bar that can accommodate 2 or 4 sling straps plus 2 chain connections
- High-performance hydraulics raise or lower individuals gradually and safely from any stationary position
- Choose either 3" or 5" casters
- Adjustable-width base

PRODUCTS

Description

13023SV	Patient Lift, New Style, Silver Vein, 1/cs
---------	--

SPECIFICATIONS

Description

Warranty	Limited Lifetime (frame), 1 Year Limited (pump)
Weight Capacity	450 lbs
Boom Height	17.5" - 72" (13023SV)
Carton Shipping Weight	84 lbs.
Casters	5" (13023SV)
Weight	72 lbs.