

## Deluxe Folding Walker, Two Button

HCPCS:



FIGURE A



- Easy push-button mechanisms may be operated by fingers, palms or side of hand
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing
- Unique U-shaped frame design creates greater clearance
- The patient is able to bring the walker closer to assist in standing because of the improved crossframe (Figure A)
- Sturdy 1" diameter, aluminum construction ensures maximum strength while remaining lightweight

### PRODUCTS

#### Description

|         |                                |
|---------|--------------------------------|
| 10200-1 | Walker, 2Btn. Adult, 1/ea      |
| 10200-4 | Deluxe Two Button Walker, 4/CS |
| 10201-4 | Deluxe Two Button Junior, 4/cs |

### SPECIFICATIONS

#### Description

|                        |                  |
|------------------------|------------------|
| Warranty               | Limited Lifetime |
| Weight Capacity        | 350 lb           |
| Carton Shipping Weight | 26   24 lb       |
| Overall Width          | 20"              |
| Weight                 | 5.5 lb           |