

Economy Folding Knee Walker

HCPCS: E0118,E0118



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E




FIGURE F



FIGURE G



FIGURE H

Keep your customers mobile while they recover from a foot or lower-leg injury with our new Folding Knee Walker. With an entry-level price point, the Knee Walker is jam-packed with all of the competitive features you're looking for, including a 350 lb. weight capacity, small footprint for optimal shelf and storage space, and retail packaging to help your customers self-educate on product features and benefits. Introduce your customers to the perfect solution for maintaining their independence and mobility!

- Stay mobile and independent during recovery
- Easily steerable, while tight turning radius allows for maneuverability in narrow spaces such as hallways and doorways
- Dual rear on-wheel braking system provides immediate stopping ability and control while navigating
- Push-button parking brake incorporated into handbrake for when you need to rest safely
- Padded leg support provides comfort, stability and is easily compatible with left or right leg
- Quickly and easily adjust the height of the kneepad or handlebar to cater to your individual needs, no tools required
- Foldable design allows for convenient storage and portability
- Strong, durable frame supports users up to 350 lb
- 8" wheels are ideal for indoor and outdoor use

PRODUCTS

Description

RTL799	Economy Folding Knee Walker, RTL, 1/EA
--------	--

SPECIFICATIONS

Description

Warranty	Limited Lifetime Lifetime Warranty on Brake Cable
Product Weight	19.5 lb
Product Dimensions	16" W x 30" L
Handle Height	31.5" - 41"
Weight Capacity	350 lb
Knee Pad Dimensions	12.5" x 7.25" x 2"
Knee Pad Height	16.75" - 23.25"
Retail Box Dimensions	20.5" (L) x 15.35" (W) x 15" (H)
Retail Box Weight	24 lb
Wheel Size	8"