

Product Catalog

LTV® Portable Ventilator





18984-001
LTV Series 1150 system

LTV Series 1150 ventilator, protective boots, AC adapter with power cord, dust cover, hand strap, circuit adapter, operator's manual, and quick reference guide



18888-001
LTV Series 1200 system

LTV Series 1200 ventilator, protective boots, AC adapter with power cord, dust cover, hand strap, circuit adapter, low flow oxygen adapter, operator's manual, and quick reference guide



22685-001
LTV2 2150 ventilator

LTV2 2150 ventilator, adapter, operator's manual, power supply, and power cord



22690-001
LTV2 2220 ventilator


LTV2 2220 ventilator, adapter, operator's manual, power supply, and power cord



18888-201
LTV Series 1200 MR Conditional ventilator system


LTV 1200 MR Conditional ventilator, protective boots, AC adapter with power cord, SprintPack™ bracket, SprintPack battery power system, dust cover, circuit adapter, low flow oxygen adapter, LTV 1200 MR Conditional floor stand, case of 5 MR-safe 15 ft circuits, and operator's manual

LTV ventilator circuits without PEEP
Single patient use (SPU) circuits, 6 ft




29657-001
Patient circuit without PEEP, 22 mm SPU (10/pkg)

1.8 m x 22 mm insp limb, wye, and 61 cm x 22 mm exp limb




29658-001
Patient circuit without PEEP, 15 mm SPU (10/pkg)

1.8 m x 15 mm insp limb, wye, and 61 cm x 15 mm exp limb



29663-001
Patient circuit without PEEP with 1 water trap, 22 mm SPU (10/pkg)

1.8 m x 22 mm insp limb with water trap, temp probe port, wye, and 61 cm x 22 mm exp limb




22685-001
Patient circuit without PEEP with 1 water trap, 15 mm SPU (10/pkg)

1.8 m x 15 mm insp limb with water trap, temp probe port, wye, and 61 cm x 15 mm exp limb


LTV circuits without PEEP (continued)

SPU heated wire circuits for Fisher & Paykel humidifiers, 5.5 ft




15090-102
Patient circuit without PEEP, 22 mm SPU single heated wire, Fisher & Paykel (10/pkg)

1.7 m x 22 mm heated insp limb, temp probe port, wye, and 61 cm x 22 mm exp limb




15090-104
Patient circuit without PEEP, 15 mm SPU single heated wire, Fisher & Paykel (10/pkg)

1.7 m x 15 mm heated insp limb, temp probe port, wye, and 61 cm x 15 mm exp limb




15091-102
Patient circuit without PEEP with elbow, 22 mm SPU dual heated wire, Fisher & Paykel (10/pkg)

1.7 m x 22 mm heated insp limb, temp probe port, wye, elbow, and 61 cm x 22 mm heated exp limb



15091-104
Patient circuit without PEEP with elbow, 15 mm SPU dual heated wire, Fisher & Paykel (10/pkg)

1.7 m x 15 mm heated insp limb, temp probe port, wye, elbow, and 61 cm x 15 mm heated exp limb



29703-001
Patient circuit without PEEP with elbow, 15 mm SPU dual heated wire, Fisher & Paykel (10/pkg)

1.7 m x 15 mm heated insp limb, temp probe port, wye, elbow, and 61 cm x 15 mm heated exp limb

SprintPack Lithium-Ion Power System



19222-201
Power system, SprintPack lithium-ion with cable protector
SprintPack power manager, cable protector, 2 SprintPack batteries, and operator's manual



19333-201
Power manager, SprintPack with cable protector
(Does not include SprintPack batteries)



19444-001
Battery, SprintPack



19096-001
Bracket, SprintPack mounting



19102-001
Transport pack for LTV and SprintPack
5 D-rings, detachable shoulder strap, hand strap and bedrail mount



11544
Cord, automotive power, 12 ft



19222-201
Mounting bracket, MR Conditional SprintPack



24555-001
SprintPack cable support bracket

External sealed lead-acid batteries



10800
Battery kit, 9-hour external

12 V/33 amp-hour battery, battery box with strap, charger, battery cable, and operator's manual



10787
Battery, 9-hour external
(12 V/33 amp-hour)



10790
Battery box with strap, external



10900
Battery kit, 3-hour external
12 V/12 amp-hour battery, padded battery bag with cable, charger, battery cable and operator's manual



11475
Battery, 3-hour external
(12 V/12 amp-hour)



10901
Cable, 3-hour external battery



11330
Bag with cable, padded battery
(for use with 3-hour battery)



10802
Cable, DC external battery
(for use with 3-hour and 9-hour battery)



10801
Charger, external battery
(for use with 3-hour and 9-hour battery)

Universal power supply (UPS)



17894-001
Battery, universal power supply

Miscellaneous power accessories



11448
AC power adapter with power cord, LTV



17527-001
AC power adapter without power cord, LTV



10854
Fuse, 32 V, 10 amp battery



11456
Electrical strip with 4 m power cord, 4-outlet U.S. hospital grade



11544
Cord, automotive power, 12 ft



10536
Power cord, 110 VAC hospital grade, 10 ft



Stands and mounting accessories



10695
Stand, LTV Series table



10611
Stand with oxygen cylinder ring, LTV Series floor

Mounts, packs and carriers




10700
Mount, bedrail




11581-G
Backpack for LTV Series ventilators with boots (gray)


Mounting accessories




10610
Bracket, female dove-tail



10979
Bracket, humidifier pole mount C-clamp



11454
Bracket, large battery mounting



17916-001
Bracket, LTM graphics monitor mounting



19096-001
Bracket, SprintPack mounting



10536
Hand strap, LTV Series ventilators (gray)



21180-001
Freedom carrier, LTV Series ventilators and SprintPack



19102-001
Transport pack for LTV Series ventilators and SprintPack
5 D-rings, detachable shoulder strap, hand strap and bedrail mount



19031-001
Hand strap, LTV Series ventilators (purple)



11501
Support arm with circuit hanger, patient circuit

Mounting accessories (continued)



11503
Hanger, clip, circuit support arm



11450
Basket, utility



11451
Crossbar, short



11452
Crossbar, long



11453
Pole, humidifer



11455
Oxygen cylinder ring set
(for use with 10611)



17341-001
Strap, transport pack shoulder

LTV Series mounting accessories sold but not shown

CAS-1152	Caster for floor stand (for use with 10611)
CAS-1168	Caster for floor stand, locking (for use with 10611)

MRI accessories



19870-001
L-bracket with screws for LTV Series
MR Conditional floor stand



19871-001
Mounting bracket, LTV Series
MR Conditional SprintPack



19873-001
Caster for LTV Series MR
Conditional floor stand




19874-001
Caster for LTV Series MR
Conditional floor stand, locking



19875-001
Oxygen cylinder ring set, LTV
Series MR Conditional floor stand

SPU MR safe circuit, 15 ft



29665-001
Patient circuit without PEEP, 22 mm SPU, 15 ft (5/pkg)
1.7 m x 15 mm heated insp limb, temp probe port, wye, elbow,
and 61 cm x 15 mm heated exp limb

Miscellaneous parts and accessories




34910
Cable, remote alarm, 50 ft
(for use with alarm 11553)



19103-001
Cable, remote alarm,
BNC - phone, 20 ft



14647-001
Cable, nurse call splitter



11550
Boot kit
(gray protective)



19175-001
Boot kit
(purple protective)



15016-001
Cap, blender



10779
Cable, nurse call, normally closed



10780
Cable, nurse call, normally open



17986-001
Adapter, electrical, Fisher
& Paykel 900 MR 900 dual
heated limb



14446
Cap tethered port
(for use with LTV Series 1150)



18098-001
Card, compact flash




10109
Cover, front panel dust



17987-001
Adapter, electrical, Fisher & Paykel
900 MR 901 single heated limb



10701
Adapter, low pressure oxygen



10702
Adapter, NIST oxygen




10678
Fan grill assembly
Fan grill, fan grill enclosure, fan filter



10258x5
Filter, LTV Series air inlet (5/pkg)

Miscellaneous parts and accessories (continued)



10789x5
Filter, LTV Series fan (5/pkg)




10699
Hose, 15 ft oxygen
(green)



11466
Test lung, 1L



10934x50
Twist-lok tie, dual (50/pkg)



10936x50
Twist-lok tie, large (50/pkg)



10696
Adapter, circuit

Ventilator shipping carton



CAT. NO.	DESCRIPTION
19818-001	Carton, shipping (replaces 17496-001)
17490-001	Sheet, divider (secures AC adapter)
17489-001	Sheet, divider (secures operator manual)
17487-001	Cap set, LTV end

Order all four parts

CAT. NO.	DESCRIPTION
19802-001	Operator's Manual, LTV 1200/1150 English
19822-001	Quick Reference Guide, LTV 1200/1150 English
18999-001	In-Service DVD, LTV 1200 English
19990-001	Caregiver Training DVD, LTV 1150 English
19960-001	Emergency Set-Up DVD, LTV 1200 English
11524	Quick Reference Guide, HME/Hospital English
21188-001	Operator's Manual, LTV 1200/1150 Spanish
21196-001	Quick Reference Guide, LTV 1200/1150 Spanish


Accessories



22770-001
Desktop charger
(English-based overlay),
batteries not included



26618-001
Battery, removable



28770-001
Battery bracket



12140-001
Patient circuit support arm



28769-001
Cord wrap bracket



29241-001
Fisher & Paykel humidifier
mounting bracket



28768-001
Basket



28767-001
Rolling floor stand
(humidifier bracket sold separately)



24425-001
Cable, nurse call, normally open
24426-001
Cable, nurse call, normally closed



28774-001
O₂ cylinder mount



30509-001
Table stand



30512-001
Bed rail mount



Need more information?

Contact Customer Service: (800) 348-9011

Terms and Conditions of Sale:

Pricing and warranties are subject to change without notice.

ZOLL Medical Corporation

269 Mill Road | Chelmsford, MA 01824 | 978-421-9655 | 800-804-4356 | zoll.com

Copyright © 2025 ZOLL Medical Corporation. All rights reserved. LTV and ZOLL are registered trademarks of ZOLL Medical Corporation in the United States and/or other countries. All other trademarks are the property of their respective owners.
MCN VP 2501 0026

ZOLL®

For subsidiary addresses and other global locations, please go to www.zoll.com/contacts.