

## Non-Tilt Overbed Tables

HCPCS:




FIGURE A



- Walnut, wood grain low-pressure laminate top
- Table top can be raised or lowered in infinite positions between 28" - 45"
- It can be raised with slightest upward pressure
- Tabletop locks securely when height adjustment handle is released
- For proper assembly, do not remove bottom screw (Figure B)
- "H" base provides security and stability
- 2" swivel casters allow for easy maneuverability
- Silver Vein (Figure A) for (Model #13067) is easy-to-maintain and less institutional looking

| PRODUCTS | Description                        |
|----------|------------------------------------|
| 13003    | Overbed Table,Dlx.Non-Tilt,1/cs    |
| 13067    | Overbed Table,SilverVein,Dlx.,1/cs |

| SPECIFICATIONS         | Description           |
|------------------------|-----------------------|
| Warranty               | Limited Lifetime      |
| Weight Capacity        | 40 lb                 |
| Base Dimensions        | 15.5" (D) x 26.5" (W) |
| Carton Shipping Weight | 24 lb                 |
| Weight                 | 21 lb                 |