

Steel Rollator with 6 Wheels, Knockdown

HCPCS: E0143 + E0156,E0143 + E0156



- Drive's Steel Rollator gives you the freedom to live a more active lifestyle without sacrificing safety or breaking the bank. Constructed of durable steel and weighing 17 lb, the rollator's reinforced frame supports users weighing up to 350 lb. For a more customized fit, easily adjust the handles ranging from 31"-36" to the desired height.
- The ergonomic hand grips allow for sustained comfort and control while loop-lock handbrakes keep the rollator safely in place. The padded seat and backrest offer a comfortable place to sit, and the flip up the seat provides easy access to a built-in storage pouch. The smooth 6" wheels allow for easy maneuverability on both indoor and outdoor terrain. When not in use, the rollator folds for compact storage and easy transport. Just lift the black strap underneath the middle of the seat, and the rollator folds flat to fit in a car trunk or closet. No tools are required, making assembly a breeze.
- Strong, durable steel frame supports up to 350 lb
- Height adjustable arms (31"-36") allow for a customized fit
- Padded seat allows user to remain comfortable during periods of extended sitting
- Padded backrest adds additional comfort and support when user is in seated position
- Convenient storage pouch under seat for easy acces to your personal items
- Ergonomic hand grip with easy to use loop-lock brakes provides comfort and control
- 6" wheels with soft grip tires ideal for indoor and outdoor use
- A compact and foldable design allows you to easily transport your rollator anywhere
- Tool-free quick and easy assembly
- Weighs just 17 lb

PRODUCTS

Description

R800KD-BL	KD Rollator, 6" 4Whl,Steel,Blue,1/ea
R800KD-GY	KD Rollator, 6" 4Whl,Steel,Gray
R800KD-RD	KD Rollator, 6" 4Whl,Steel,Red,1/ea

SPECIFICATIONS

Description

Product Weight	17 lb
Seat Dimensions	12" (W) x 12" (D)
Handle Height	31" - 36"
Weight Capacity	350 lb
Seat To Floor Height	20"
Carton Shipping Weight	19.84 lb
Dimensions	24.8" (L) x 22.25" (W) x 31.5" (H) Basket: 7" (L) x 12" (W) x 5" (H)
Wheel Size	6"
Width Between Arms	18.5"
Width inside back legs	21.6"