

Transfer Tub Bench

HCPCS: E0247,E0247



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Tool-free assembly of back, legs and arm
- A-frame construction provides stability
- Durable blow-molded plastic bench and backrest
- Height adjusts in ½" increments with unique "Dual Column" extension legs
- Reversible to accommodate any bathroom
- Pinch-free lever allows for push pins to be depressed without pinching fingers
- Made using recycled materials

PRODUCTS	Description
12011KD-1	Transfer Bench KD,Tool Free, 1/cs
12011KD-2	Transfer Bench KD Tool Free,2/cs
SPECIFICATIONS	Description
Warranty	Limited Lifetime
Seat Dimensions	26" (W) x 18.5" - 19.5" (D)
Weight Capacity	400 lb
Carton Shipping Weight	13.6 lb
HCPCS	E0247
Height Adjustment	17.5" - 22.5"
Weight	13.6 lb
ACCESSORIES	Description
12011SS	SuctionCups,Sm.12011,1pr/bx, BX
12011SL	SuctionCups,Lg.12011,1pr/bx, BX