

Steel Transport Chair

HCPCS: E1038




FIGURE A



FIGURE B

- Weighs 26 lb
- Attractive silver vein finish
- Available in 17" or 19" seat widths
- Back folds down for storage and transport
- Swing-away footrests come standard
- Padded armrests provide additional comfort
- 8" casters in front and rear
- Comes with attractive nylon upholstery that is durable and easy-to-clean
- Standard rear wheel locks
- Comes with upholstery handles to fold chair
- Convenient cup holder can easily attach to chair Sold separately (Model# STDS1040S)

PRODUCTS

Description

TR37E-SV	Trans Chair Econ 17in Slv Vein 1/cs
TR39E-SV	TransChair,Econ 19in SlvVein,1/cs

SPECIFICATIONS

Description

Warranty	Limited Lifetime
Weight Capacity	250 lb
Seat To Floor Height	19"
Carton Shipping Weight	31 lb
Weight	26 lb
Back Height	18"
CA Prop 65	[WARNING: This product can expose you to chemicals that are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to www.P65Warnings.ca.gov .]

ACCESSORIES

Description

STDS1040S	Cup Holder, Swivel,1BX,6/CS, BX
-----------	---------------------------------